

IPHONE 7.4 UPGRADE INSTRUCTIONS 5/10/12

To delete the current version press and hold the PatientKeeper icon until the screen shakes.



After the screen shakes the icons will begin to wave side to side and have an "x" next to them.



Click the X next to the PatientKeeper icon (PK53154R).




When prompted to delete PK53154R click Delete.




The main screen will be displayed with the wavy icons and the delete option. Press the menu button on the phone to regain normal operation.




Open the App store and do a search for PK7426R



Select the PK7426R app



Click Install.




Once the app is installed launch the app by clicking PK7426R.




Enter and save the host name (A list of host names can be found in the Portal Device Connectivity Guide).

South – dedv.mobile.medicity.net

North – mwdv.mobile.medicity.net



Log in to PK.



The screenshot shows the mobile app interface for PatientKeeper. At the top, the status bar displays 'AT&T', signal strength, Wi-Fi, the time '10:15 AM', and battery level. Below the status bar is a navigation bar with a 'Server Info' button on the left and the title 'Login' in the center. The main content area features the 'patientkeeper' logo, followed by two text input fields labeled 'username' and 'password'. A blue 'Login' button is positioned below the password field.

User will be required to provision their device on the 1st login. After provisioning the device will sync and the patient list will be displayed.

